

MTL.12.100&200
Inverted Metallurgical Microscope



Metallurgical Microscope Outfits

Items	Specifications	MTL.12. 100	MTL.12. 200
Eyepieces	WF10×/Φ18	●	●
	WF10×/(Reticule, 0.1mm)	○	●
Metallurgical LWD Infinity Plan Objectives	LPL 5× / 0.13 W.D=16.04mm	●	●
	LPL 10× / 0.25 W.D=18.48mm	●	●
	LPL 20× / 0.40 W.D=8.35mm	○	○
	LPL 50× / 0.70(S) W.D=1.95mm	●	●
	LPL 80× / 0.80(S) W.D.= 0.85mm	○	○
	LPL 100× / 0.85(Dry) W.D=0.3mm	○	○
Binocular Head	Inclined 45°	-	●
Monocular Head	Inclined 45°	●	-
Nosepiece	Quadruple	●	●
Mechanical Stage	Size: 180×180mm Travel Range: 30X30mm	●	●
Stage Plate	1 (Φ10)	●	●
	2 (Φ20)	●	●
	3 (Φ40)	○	○
Set of Stage Plate		●	●
Main Body		●	●
6V/20W Halogen Lamp	Input voltage: 100V~240V	○	○
3W-LED Lamp		●	●
Filters	Blue, Green, Grey	●	●
Micrometer	Division 0.01mm	○	○

Note: "●" In Table Is Standard outfits, "○" Is Optional Accessories.

Packing Size: 437mm×268mm×522mm Gross Weight: 8 kgs Net Weight: 7 kgs

Metallurgical Microscope Size (mm)

